



THUASNE®


LombaStab High

fr	Ceinture de soutien dorso-lombaire - "Serrage facile et précis"*	5
en	Dorso-lumbar support belt - "Easy, precise tightening"*	6
de	LWS-Stabilisierungsorthese - „einfaches und präzises Justieren“*	7
nl	Dorso-lumbale rugbrace - "Eenvoudige en nauwkeurige aanspanning"*	8
it	Cintura di sostegno dorso-lombare - "Serraggio facile e preciso"*	10
es	Faja de sujeción dorsolumbar - "Fijación fácil y precisa"*	11
pt	Cinta de suporte dorso-lombar - "Aperto fácil e preciso"*	12
da	Rygstøttebælte - "Nem og præcis stramning"*	14
fi	Selkä- ja lantiotukivyo - "Helpostija tarkasti kiristettävä"*	15
sv	Ryggstöd för ländryggen - "enkel och precis åtdragnings"*	16
el	Zώνη οσφυϊκής υποστήριξης - «Εύκολη και ακριβής σύσφιξη»	17
cs	Dorso-lumbální podpůrný pásek - „Jednoduché a přesné stažení“*	19
pl	Orteza lędźwiowo-krzyżowa - „Łatwe i precyzyjne mocowanie“*	20
lv	Krušu-jostasvietas atbalsta josta - „Viegla un precīza savilkšana“*	21
lt	Palaikomasis krūtinės ir juosmens diržas - „Lengvai ir tiksliai suveržiamas“*	22
et	Selja-ja nimme piirkonna tugivõö - „Lihitne ja täpne pingutus“*	24
sl	Hrbtno-ledvena opornica - „Enostavno in natančno vpetje“*	25
sk	Podporný chrabtovo-driekový - „jednoduché a presné dotiahnutie“*	26
hu	Háti-ágyékű fűző - „Egyszerű és pontos rögzítés“*	27
bg	Дорзолумбален колан - „лесно и точно затягане“*	28
ro	Centură de susținere dorso-lombară - „Strângere usoară și precisă“*	30
ru	Опорный поясничный - «Легко и просто затянуть»*	31
hr	Lumbalno-leđni potporni pojasi - „Jednostavno i precizno stezanje“*	33
zh	调节简易准确的胸腰部支撑带*	34
ar	حزام لدعم الصدر وأسفل الظهر - "شد سهل ودقيق".....	37

Instructions d'utilisation - Instructions for use - Gebrauchsweisung - Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso - Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet - Råd vid användning - Oδηγίες χρήσης - Pokyny k používání - Szczegóły użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo - Pokyny na používanie - Használati utasítás - Инструкции за употреба - Instrucțiuni de utilizare - Инструкция по эксплуатации - Upute za uporabu - **使用指南** - ارشادات الاستعمال

	1	60 - 80 cm
	2	75 - 95 cm
	3	90 - 110 cm
	4	105 - 125 cm
	5	120 - 140 cm

fr	Tour de taille	Hauteurs : 35 cm - 40 cm
en	Waist circumference	Heights: 35 cm - 40 cm
de	Taillenumfang	Höhen: 35 cm - 40 cm
nl	Tailleomvang	Hoogtes: 35 cm - 40 cm
it	Circonferenza vita	Altezze: 35 cm - 40 cm
es	Contorno de cintura	Alturas: 35 cm - 40 cm
pt	Perímetro de cintura	Alturas: 35 cm - 40 cm
da	Taljemål	Højde: 35 cm - 40 cm
fi	Vyötärönympärys	Korkeudet: 35 cm - 40 cm
sv	Midjemått	Höjd: 35 cm - 40 cm
el	Περίμετρος μέσης	Ύψη: 35 εκ - 40 εκ.
cs	Obvod pasu	Výšky: 35 cm - 40 cm
pl	Obwód w talii	Wysokości: 35 cm - 40 cm
lv	Vidukļa apkārtmērs	Augstumi: 35 cm - 40 cm
lt	Juosmens apimtis	Aukštis: 35 cm - 40 cm
et	Võõrumbermõõt	Kõrgused: 35 cm - 40 cm
sl	Obseg pasu	Višina: 35 cm - 40 cm
sk	Obvod pása	Výšky: 35 cm - 40 cm
hu	Derékbőség	Magasság: 35 cm - 40 cm
bg	Обиколка на талията	Височини: 35 см - 40 см
ro	Circumferința taliei	Înălțimi: 35 cm - 40 cm
ru	Обхват талии	Высота: 35 см - 40 см
hr	Opseg struka	Visine: 35 cm - 40 cm
zh	腰围	高度 : 35 cm - 40 cm
ar	محيط الخصر	الارتفاعات: 35 سم - 40 سم

fr	Stabilisation	Éducation posturale	Effet proprioceptif
en	Stabilization	Postural education	Proprioceptive effect
de	Stabilisierung	Erlernen der Körperhaltung	Propriozeptive Wirkung
nl	Stabilisatie	Voorlichting over houding	Proprioceptieve werking
it	Stabilizzazione	Educazione posturale	Effetto propriolettivo
es	Estabilización	Educación postural	Efecto proprioceptivo
pt	Estabilização	Educação postural	Efeito propriojetivo
da	Stabilisering	Postural genoptræning	Proprioceptisk effekt
fi	Stabilointi	Asentokuntoutus	Liikeaistiin kohdistuvaa vaikutusta
sv	Stabilisering	Postural tränning	Proprioceptiv effekt
el	Σταθεροποίηση	Σωστή στάση σώματος	Αποτέλεσμα ιδιοδεκτικότητας
cs	Stabilizace	Správné držení těla	Proprioceptívní účinek
pl	Stabilizacja	Korekcja postawy	Propriocepcaja
lv	Stabilizēšana	Stājas korekcija	Proprioceptīvs efekts
lt	Stabilizavimas	Laikysenos taisymas	Propriocepcinis poveikis
et	Stabiliseerimine	Rühitreening	Propriiotseptiivne toime
sl	Stabiliziranje	Učenje telesne drže	Proprioceptivni učinek
sk	Stabilizácia	Cvičenie posturálnych svalov	Proprioceptívny účinok
hu	Stabilizálás	Tartásjavító oktatás	Proprioceptív hatás
bg	Стабилизация	Приучаване към правилна стойка	Проприоцептивен ефект
ro	Stabilizare	Educație posturală	Efect proprioceptiv
ru	Стабилизация	Исправление осанки	Проприоцептивный эффект
hr	Stabilizacija	Učenje o držanju	Učinak proprioceptivnog vježbanja
zh	稳定	姿势联系	本体感
ar	الاستقرار	تدريب وضعی	مفعول استقرار الحس العميق

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Conserver cette notice.

en

DORSO-LUMBAR SUPPORT BELT - "EASY, PRECISE TIGHTENING"*

Description/Destination

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

Composition

Elastic fabric. ①

Handles. ②

Back panel with Coolmax lining. ③

Abdominal panels with finger-loops. ④

Laces.

Textile components: polyamide - polyester - elastane - viscose - polyurethane.

Non-textile components: brushed tempered steel - polypropylene - SEBS - polyamide - polyurethane - polyoxymethylene.

Properties/Mode of action

The device relieves dorso-lumbar pain by providing support to the lumbar spine: compartmental lumbar

unloading is therefore achieved via increased intra-abdominal pressure coupled with a 3-pressure-point system.

This device uses back pads, implementing a pulley-like system to distribute pressure.

Indications

Non-specific low back pain (acute, subacute and chronic phases).

Spinal stenosis.

Spondylolisthesis/spondyloarthritis.

Stable vertebral fracture (osteoporotic).

Post-operative immobilisation.

Osteochondrosis/Scheuermann's disease.

Spondylolistesis.

Contraindications

Do not use in the event of known allergy to any of the components.

Do not use on pregnant women.

Do not use in the event of an unstable vertebral fracture.

Do not use in the event of hiatus or abdominal hernia, or of severe liver disease.

Do not use in the event of bone cancer with spinal metastases.

Do not use in the event of circulatory, pulmonary or cardiovascular problems in patients for whom an increase in arterial blood pressure is not recommended.

Precautions

Verify the product's integrity before every use.

Wear the device above thin clothing.

Do not use the device if it is damaged.

Choose the appropriate size to fit the patient, referring to the size chart.

Strictly comply with your healthcare professional's prescription and recommendations for use.

This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

Do not wear the device while sleeping.

In the event of discomfort, significant hindrance, pain, abnormal sensations, remove the device and consult a healthcare professional.

The use of a medical device by a child should be done under the supervision of an adult or a healthcare professional.

Measure the waist circumference of the patient at the mid-distance between the iliac crests and floating ribs.

It is recommended that a healthcare professional supervises the first application.

For hygiene, security and performance reasons, do not re-use the product for another patient.

For safety reasons, do not wear the belt close to machines (risk of becoming caught).

Do not apply the product in direct contact with skin (risk of pinching).

Do not wear the device in a medical imaging machine.
Do not wear the device when driving a vehicle.

Undesirable side-effect

This device can cause skin reactions (redness, itching, irritations...).

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

Instructions for use/Application

It is recommended that a healthcare professional adjust the stays to the shape of the patient's back before the belt is worn for the first time.

Check that the two tightening handles ❷ are fully retracted.

Open the belt and lay it out completely flat so that you can see the handles. ❸

If necessary, unhook them and separate the two back parts by pulling the ends of the belt, then fasten the handles again.

Position the belt around the waist, using the finger loops located on the front panels. ④

The tightening back pads must be positioned outwards (i.e.: not against the skin).

The lacing section must be centered over the spine: the pads should be positioned on either side of the spine. ⑤

The bottom of the belt must stop at the top of the crease between the buttocks. ⑥

Thusne® logo indicates which way up the belt should be put on.

Fasten the belt: right side over left side, without pulling too tight or applying any particular pressure. ⑦ ⑧

Adjust the tightness using the two handles: pull them forwards then fold them flat against the front of the belt, where they attach like a self-fastener. ⑨ ⑩

These two handles can be used to adjust the level of compression as required throughout the day.

Care

Close the self-fastening tabs before washing. Machine washable at 30°C (delicate programme). Do not machine wash the device more than 10 times. Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Do not tumble-dry. Do not iron. Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

Storage

Store at room temperature, preferably in the original packaging.

Disposal

Dispose of in accordance with local regulations.

*Study conducted in-house on a panel of 13 people, July 2019.

This medical device is a regulated health product with CE marking.

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©Thuasne - 2054201 (2024-04)
THUASNE UK Ltd
Thuasne SAS – SIREN/RCS Nanterre 542 091 186