THUASNE®

# LombaStab Dorso

fr	Ceinture de soutien thoraco-lombaire - "Serrage facile et précis"*	5
en	Thoracolumbarsupport belt - "Easy, precise tightening"*	6
de	BWS-/LWS Stabilisierungsorthese - „Einfaches und präzises Justieren“*	8
nl	Thoraco-lumbale rugbrace - "Eenvoudige en nauwkeurige aanspanning"*	9
it	Cintura di sostegno toraco-lombare - "Serraggio facile e preciso"*	11
es	Faja de sujeción toracolumbar - "Fijación fácil y precisa"*	13
pt	Cinta de suporte toracolombar - "Aperto fácil e preciso"*	14
da	Lændestøttebælte med thorakal støtte - "Nem og præcis stramning"*	16
fi	Rinta- ja lannerangan tukivyo - "Helposti ja tarkasti kiristettävä"*	17
sv	Ryggstöd för bröstryggen - "Enkel och precis åtdragning"*	19
el	Zώνη θωρακο-οσφυϊκής υποστήριξης - «Εύκολη και ακριβής σύσφιξη»*	20
cs	Pás pro podporu hrudní a bederní páteře - „Jednoduché a přesné stažení“*	22
pl	Orteza piersiowo-lędźwiowo-krzyżowa - „Łatwe i precyzyjne mocowanie“*	24
lv	Krūšu-jostasvietas atbalsta josta - „Viegla un precīza savilkšana“*	25
lt	Palaikomasis krūtinės ir juosmens diržas - „Lengvai ir tiksliai suveržiamas“*	27
et	Nimme piirkonna tugikorsett - „Lihtne ja täpne pingutus“*	28
sl	Oporni prsno-križni pas - „Enostavno in natančno vpetje“*	30
sk	Hrudno-driekový podporný pás - „Jednoduché a presné dotiahnutie“*	31
hu	Háti-ágyéki fűző - „Egyszerű és pontos rögzítés“*	33
bg	Торако-лумбален колан - „лесно и точно затягане“*	34
ro	Centură de susținere toraco-lombară - „Strângere ușoară și precisă“*	36
ru	Грудопоясничный ортез с системой простой и быстрой фиксации*	37
hr	Torakalno-lumbalni potporni pojaz - „Jednostavno i precizno stezanje“*	39
zh	调节简易准确的胸腰部支撑带*	41
ar	حزام لدعم الصدر وأسفل الظهر - "شد سهل ودقيق".....	43

Instructions d'utilisation - Instructions for use - Gebrauchsanweisung - Gebruiksaanwijzingen - Istruzioni per l'uso - Instrucciones de uso - Instruções de utilização - Anvisninger vedrørende brug - Käyttöohjeet - Råd vid användning - Одынгес хронас - Pokyny k používání - Szczegóły użytkowania - Lietošanas instrukcija - Naudojimo instrukcijos - Kasutusjuhised - Navodila za uporabo - Pokyny na používanie - Használati utasítás - Инструкции за употреба -

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Conserver cette notice.

**en**

## **THORACOLUMBARSUPPORT BELT - "EASY, PRECISE TIGHTENING"\***

#### **Description/Destination**

The device is intended only for the treatment of the indications listed and for patients whose measurements correspond to the sizing table.

#### **Composition**

Elastic fabric. ①

Adjustable back frame. ②

Straps with hooks. ③

Sternum support. ④

2 abdominal panels with 2 finger-loops. ⑤

Handles. ⑥

Back panel with Coolmax lining.

Laces.

Textile components: polyamide - polyester - elastane - viscose - polyoxymethylene - polyethylene - elastodiene.

Non-textile components: aluminium - polypropylene - SEBS - polyamide - polyurethane.

#### **Properties/Mode of action**

The device relieves lower back pain by supporting the spine through an increase in intra-abdominal pressure. This reduces the load on the intervertebral discs.

This device uses back pads, implementing a pulley-like system to distribute pressure.

The device maintains spinal hyperextension thanks to the adjustable straps, comfortable metal back frame and removable sternum support.

#### **Indications**

Stable vertebral fracture (traumatic and osteoporotic). Scheuermann's disease.

Static disorders (kyphosis, lordosis).

Non-specific low back pain (acute, subacute and chronic phases).

#### **Contraindications**

Do not use in the event of known allergy to any of the components.

Do not use on pregnant women.

Do not use in the event of an unstable vertebral fracture.

Do not use in the event of hiatus or abdominal hernia, or of severe liver disease.

Do not use in the event of bone cancer with spinal metastases.

Do not use in the event of circulatory, pulmonary or cardiovascular problems in patients for whom an increase in arterial blood pressure is not recommended.

#### **Precautions**

If necessary, the healthcare professional who prescribed or delivered the device can cut the frame pouch and straps to the desired length with scissors.

Verify the product's integrity before every use.  
Wear the device above thin clothing.  
Do not use the device if it is damaged.  
Choose the appropriate size to fit the patient, referring to the size chart.  
Strictly comply with your healthcare professional's prescription and recommendations for use.  
This product is intended for the treatment of a given condition. Its duration of use is limited to this treatment only.

Do not wear the device while sleeping.

In the event of discomfort, significant hindrance, pain, abnormal sensations, remove the device and consult a healthcare professional.

The use of a medical device by a child should be done under the supervision of an adult or a healthcare professional.

Measure the waist circumference of the patient at the mid-distance between the iliac crests and floating ribs. It is recommended that a healthcare professional supervises the first application.

For hygiene, security and performance reasons, do not re-use the product for another patient.

Do not apply the product in direct contact with skin (risk of pinching).

For safety reasons, do not wear the belt close to machines (risk of becoming caught).

Do not wear the device in a medical imaging machine.  
Do not wear the device when driving a vehicle.

#### **Undesirable side-effect**

This device can cause skin reactions (redness, itching, irritations...).

Any serious incidents occurring related to the device should be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is resident.

#### **Instructions for use/Application**

It is recommended that a healthcare professional adjust back frame to the shape of the patient's back before first use.

#### **Back frame adjustment**

Measure the height of the patient's spine such that the lower edge of the back frame is level with the coccyx and the upper edge is level with thoracic vertebra T1 .

Loosen the four screws slightly and adjust the frame height before tightening the screws.  

Without using tools or heat source, adjust the frame to the curve of the patient's back.

Insert the fitted frame into the back pouch.

Shorten the pouch to the defined frame length by cutting it with scissors. 

#### **Strap adjustment**

Ask the patient to tighten the straps until they are flat against his/ her body, such that they straighten the back slightly and stabilise the spine.

Shorten the straps used to tighten the back section using scissors  .

Fit the hooks back onto each strap.

#### **Fitting the device**

Open the belt and lay it out completely flat so that you can see the handles .

Check that the two tightening handles are fully retracted

If necessary, unhook them and separate the two back parts by pulling the ends of the belt, then fasten the handles again.

Put on the brace in the same manner as a backpack .

Position the belt around the waist, placing your hands in the finger loops  located on the abdominal panels .

The tightening back pads must be positioned outwards (i.e. : not against the skin).

The lacing section must be centered over the spine: the pads should be positioned on either side of the spine.

The bottom of the belt must stop at the top of the crease between the buttocks.

The Thuisne® logo indicates which way up the belt should be put on.

Fasten the belt: right side over left side, without pulling too tight or applying any particular pressure .

Adjust the tightness using the two handles : pull them forwards then fold them flat against the front of the belt, where they attach like a self-fastener ( .

Position the sternum support  and adjust the tightness; it should be positioned over the sternum, above the chest, without compressing the thorax .

Adjust strap tension  and place the hooks on the self-fastening surface ( 

Avoid placing the hooks against the tightening system thread as they can cause damage.

#### **Care**

Remove the back frame before washing. Close the self-fastening tabs before washing. Machine washable at 30°C. Do not machine wash the device more than 10 times. Do not dry clean. Do not use detergents, fabric softeners or aggressive products (products containing chlorine, etc.). Do not tumble-dry. Do not iron. Squeeze out excess water. Dry away from any direct heat source (radiator, sun, etc.). Dry flat.

#### **Storage**

Store at room temperature, preferably in the original packaging.

#### **Disposal**

Dispose of in accordance with local regulations.

\*Study conducted in-house on a panel of 13 people, July 2019.

This medical device is a regulated health product with CE marking.

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Den Beutel abhängig von der festgelegten Gestelllänge mit einer Schere kürzen. ©
**THUASNE UK Ltd**
Thuasne SAS – SIREN/RCS Nanterre 542 091 186